NOHRD Seated Row

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The screenshot shows a training interface with a target area at the top. The target area has a 'Ziel' (Target) value of 15 and 'Wiederholungen' (Repetitions) of 0. A black dot represents the cursor. Below the target area, there is a brown bar with the text: 'Führen Sie 3 Wiederholungen durch, um Ihren Bewegungsbereich zu kalibrieren.' Below this bar, there are three small circles. At the bottom, there is a table with four rows: 'Verstrichene Zeit' (00:00), 'Bewegungsbereich' (-), 'Exzentrische Zeit' (-), and 'Konzentrische Zeit' (-). A red 'STOPPEN' button is at the bottom.

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GENERAL INFORMATION

DE The NOHRD Seated Row is designed for versatile back training using
EN cable-based exercise. The equipment is fitted with a newly developed lower-body support system for added stabilization during training.

Intended Use

The equipment complies with Class S + I in accordance with DIN EN ISO 20957-1 and is intended for use in leisure sports facilities and therapeutic institutions where access and supervision are regulated by the owner or operator. The level of supervision required depends on the user's reliability, age, experience, and overall condition.

Individual Adjustment

No pre-setting of the equipment is required prior to training. Training resistance is adjusted individually by selecting the desired weight plates.

Always ensure that the selector pin is fully inserted so that the magnetic locking mechanism engages securely.

SEATED ROW

Starting Position

Grip both handles and sit down on the seat.

Position the knees in front of the knee pad and assume an upright seated posture.

Movement

From the starting position, pull the slightly bent arms back toward the body, then extend them again in a controlled manner.

Keep the thumbs pointing upward throughout the movement.

Target Muscles

Latissimus dorsi, shoulder girdle muscles, and elbow flexors.

DE
EN


INSTALLATION

DE
EN

Position the equipment on a stable and level surface.

Any floor unevenness can be compensated using the adjustable leveling unit located beneath the base frame.

Ensure that sufficient safety clearance is maintained from surrounding objects, such as other equipment or walls.

Technical Specifications

Dimensions: L 163 / W 95 / H 170 cm

Footprint: approx. 1,
m²

Equipment weight: approx. 151 kg

Maximum training load: approx. 75 kg

Maximum user weight 135 kg

Approved for use by 1 person.

WARNING NOTICE

Training Area

Training is permitted exclusively on the equipment. During all exercises, your center of gravity must remain within the designated equipment area. Training in a lateral position next to the equipment is strictly prohibited.

DE
EN

Unsupervised Children


Unsupervised children must be kept away from the equipment at all times.

Environmental Conditions

Transport, storage, and operation of the NOHRD Seated Row in damp or humid environments are not permitted. During operation, the ambient temperature must be between 5–40 °C.

Removal of the Wooden Covers

To perform maintenance work, the rear wooden panel must be removed. Loosen the securing screws on the side frame. The rear panel can then be slid upward and lifted off.

Monitor Operation (optional accessory)

DE
EN

The monitor activates automatically when the equipment is connected to the power supply. The equipment does not feature a separate On/Off switch.

Wi-Fi Connection

The monitor software is continuously updated and improved.

When restarting the monitor, the latest software version will be installed automatically, provided that the monitor is connected to your Wi-Fi network.

To connect your monitor to Wi-Fi, follow the steps below:

Once the monitor has completed startup, draw a large "M" on the touchscreen using your finger (see image).


A window will now open, prompting you to enter a code.

The code for your Seated Row is:
1708

IMPORTANT:

Do not discard this instruction manual. You can refer to it at any time to retrieve the code.

Alternatively, you may contact our customer service team to obtain the code.


You are now in the settings menu.

Under “Network & Internet,” select and configure your Wi-Fi network.

Once the Wi-Fi connection has been successfully established, swipe the settings menu upward from the bottom of the screen to return to the home screen.

IMPORTANT:

Never perform a “Factory Reset.” This will erase the program and can only be reinstalled at the factory.

DE
EN


The home screen provides options for language selection as well as two training modes:

a) Quick Start

b) Exercises


This mode displays exercises that can be performed using the equipment.

Select the desired training mode.

If you select "Exercises," you will be directed to a selection menu displaying multiple exercises.

Select the desired exercise.

After answering a few questions regarding your gender and your self-assessed physical condition, the equipment will suggest a recommended training weight.

CAUTION:

This recommendation is based on scientific principles but may not be suitable for every individual.

Before you begin training, instruction by a certified trainer is strongly recommended.

Press START to begin.

←

READY TO START?
Confirm your selected settings before starting the session

Exercise	Level 1 - Pull-down
Sex	Male
Experience	Advanced
Suggested weight	8 plates (36 kg)

Your equipment weight is personalised based on your sex and experience level.

START

You will now be prompted to perform three repetitions of the selected exercise in order to determine your individual range of motion together with your trainer.

The system calculates the average of these three repetitions.

target 15 Repetitions 0

Perform 3 repetitions to calibrate your range of motion.

elapsed time	movement range
00:00	-
eccentric time	concentric time
-	-

STOP


DE
EN

A cursor on the right side of the screen will move in sync with your motion.


Your objective is to reach the highlighted area at both the top and bottom of the movement bar.

Once this target is achieved, one repetition is counted.

After each completed range of motion, the accuracy of the movement (based on the average of the first three repetitions) is displayed as a percentage.

Upon completion of the exercise, a summary screen will be displayed.

WARNING

Incorrect or excessive training may result in serious injury or health risks. The recommended exercises for the NOHRD Seated Row are intended exclusively for healthy individuals. If in doubt, consult your physician to determine whether you are fit to train with this equipment. Before beginning any training session, ensure that the equipment is in proper working condition. The equipment may only be operated with the factory-installed additional weights.

Caution

Ensure proper movement execution during this exercise. As a general rule, use light weights.

CARE AND MAINTENANCE

The safety level of the equipment can only be maintained if it is inspected regularly for damage and wear. This includes adjustment mechanisms, locking screws, and all moving parts. Other components subject to wear include load-bearing elements such as cables and carabiners. Defective parts must be replaced immediately with original replacement parts, or the equipment must be taken out of service until repairs have been completed.

Frame components and upholstery should be cleaned using water only. Individuals with open wounds must not train on the equipment, as it may only be cleaned using water-soluble disinfectants.

WATERROWER | NOHRD GmbH • Otto-Hahn-Straße 75 • 48529 Nordhorn • Germany

Tel. +49 (0) 59 21 - 17 98 400 • Fax +49 (0) 59 21 - 17 98 411

info@nohrd.de • www.nohrd.com